

OctaFlow™ II

MultiValve, MultiFunction Superfusion/Perfusion System

OctaFlow™ II: ALA's most versatile Drug Application Device for electrophysiology and imaging research.



The **OctaFlow™ II** offers the following user benefits:

- Modular banks of reservoirs for up to 32 samples per experiment
- Choice of valves for rapid exchange or minimal maintenance
- Valve-control of solution flow instead of motorized manipulation of barreled pipettes, minimizing sample losses following solution exchange
- Keyboard solution selection, for exchange “on the fly”
- Trigger in/out for synchronization with 3rd party data acquisition systems
- Preprogrammed solution selection from GUI or manual
- Choice of reservoir sizes (1ml, 5ml, 10ml, or 60ml; 5ml standard). (glass options available)
- Rapid flush mode for easy cleaning of system
- Pressurization of fluids for optimal output
- Programmable analog output reports valve identification, step number or pressure to data acquisition system
- USB 2.0 compatible for fast PC communication
- Windows 10 compatible
- Millisecond solution exchange time




OctaFlow™ II, Software



The **OctaFlow™ II** is a fully software-controlled system. The software provides five areas of functionality detailed below.

- 1 Valve and reservoir configurations can be saved for future experiments. Customizations include unique pressures and labels for individual reservoirs as well as unique voltages for individual valves.
- 2 Perfusion sequence setup includes parameters such as valve open times, duration, flush time, delay time, and output trigger in a flexible, easy-to-use spreadsheet format.
- 3 Program run control allows for initiation and termination of protocols while monitoring progress and configuration settings.
- 4 Manual control can be achieved through individual activation. Hot keys are integrated for on-the-fly control without protocol setup. Manual sequences can be saved for use as future macros.
- 5 Additional tools, available for greater scope of performance, initiate automated cleaning sequences, select global software settings, and create log files.

OctaFlow™ II software offers "on the fly" control of pressure and times settings, switches from single- to dual-valve mode, analog output control, input/output triggers, and lots more!



Sample pressurized flow rate:
1ml in 9 minutes at 10 PSI/ 70kPa through QMM Micro-manifold® with 100µm ID tubes

Ordering information key code:

OctaFlowXYZ:

XX = # of valves: specify 8, 16, 24, or 32 valves;

Y = type of valve: order **S** for custom Lee solenoid valve or **P** for 3-way pinch valve;

Z = # of **QMM MicroManifolds®** included - Order **1**, **2**, **3** or **4** where:

1 = **QMM** type is same as the number of valves ordered - **1 QMM** included with 1 flush valve

2 = **QMM-16** is included as base and the other is **QMM-8** for 24 valves or **QMM-16** for 32 valves with 2 flush valves

3 = **QMM-8** is included as base and then 2 more **QMM-8**'s for 24 valves or 1 more **QMM-8** and 1 more **QMM-16** for 32 valves with 3 flush valves.

4 = 4 x **QMM-8**'s are included with 4 flush valves.

Examples of ordering code:

OctaFlow32P/4: 32 channel pinch valve system with 4x **QMM-8 MicroManifolds®** and 4 flush valves.

OctaFlow24S/2: 24 channel solenoid valve system with 1x **QMM-16**, 1 x **QMM-8 MicroManifolds®**, and 2 flush valves.

* **OctaFlow 32P/1 & 32S/1** are not available. All other 32 channels are available.

OctaFlow™ Specifications*:

Max. # of perfusion valves: 32	Max. suction developed via internal Venturi pump: 199 mmHg nominal
Max. # of valves/Bank: 8 valves/bank, up to 4 banks	Typical speed of pressure rise: 520mmHg/sec
Max. current output per valve: 100mA (Lee Solenoid) 350mA (Pinch)	Sample flow rate: 1ml in 9 min. @ 520 mmHg w/standard QMM
Max. Voltage Per Valve: 18 volts	Standard QMM: 8 tubes @ 100µm ID, 9 th tube @ 200µm ID, 100µm ID tip
Max. ontime per valve: no limit Min. ontime per valve: 2ms @ 12V & 1ms @ 18V	USB 2.0: Support Software is compatible with Windows 7 & 10 (32-& 64-bit)
Max. input pressure: 3970mmHg - other units selectable	Power requirements: 120/220V
Max. pressure applied: ~800mmHg	OctaFlow interface: 19"x3.5"x10" rack mountable, 10 lbs/4.5kg
Max. No. of Sequence Steps: 254	Programmable Voltage Range: 1-18volts/valve
* specifications are subject to change without prior notifications	



Scientific
Instruments

Web: www.alascience.com
E-Mail: sales@alascience.com