

## Technical Specifications

### EPC 7 PLUS Main Unit:

#### Dimensions Main Unit (D x W x H):

(31.1 x 48.3 x 14.5) cm,  
(12.3 x 19.0 x 5.7) inch

#### Weight:

7.6 kg / 16.5 lbs

Mounts in a 19" rack.

Operates on standard 115 V / 230 V.

### EPC 7 PLUS Headstage:

#### Dimensions (D x W x H):

(70 x 40 x 19 ) mm  
(2.75 x 1.57 x 0.75 ) inch

#### Current measuring resistors:

50 G (high range)

500 M (medium range)

#### Largest measurable currents:

200 pA (high range)

20 nA (medium range)

#### Noise measured with open input:

(8-pole Bessel filter, high range)

up to 1 KHz: < 30 fA RMS

up to 3 KHz: < 85 fA RMS

up to 10 KHz: < 350 fA RMS

### Current Monitor Signal:

Gain: 0.5 to 1000 mV/pA, switch-selected

#### Maximum bandwidth:

100 KHz (medium range)

60 KHz (high range)

### Filters:

The Current Monitor signal of the EPC 7 PLUS is filtered by a 7-pole Bessel filter from 0.1 to 30 kHz in 11 steps. In the additional 'FULL' setting, a 3-pole Bessel filter is used and the current signal is provided at full bandwidth of the EPC 7 PLUS.

### Capacitance Compensation:

Manual adjustment of the fast and slow capacitance cancellation.

*C-fast:* 0 to 10 pF, 0.5 to 5  $\mu$ s tau

*C-slow:* 0.2 to 10 pF, calibrated

2 to 100 pF, calibrated

10 to 1000 pF

#### Series conductance adjustment:

0.01 to 1  $\mu$ S

### Series Resistance Compensation:

Automatically determined from transient cancellation controls.

*Fractional compensation:* 0 – 90 %

### Pipette Offset:

Manual adjustment within a  $\pm$  200 mV range.

### Holding Potential:

Hardware controlled holding potential with a  $\pm$  500 mV range.

### Current Clamp CC Command Resolution:

1 pA/mV input; up to  $\pm$ 1 nA

### Noise Monitor Facility:

LCD displays RMS current at gain setting of 100 mV/pA in units of fA for checking background noise. The bandwidth is set by the filter frequency.

## Related Products

### LIH 1600

16-bit multi-channel data acquisition system.

### PATCHMASTER

Multi-channel patch clamp data acquisition and analysis software for Windows and Mac OS.

### FITMASTER

Data analysis software supporting PATCHMASTER data structure for Windows and Mac OS.

### PULSE/PULSEFIT

Patch clamp data acquisition and analysis software for Windows and Mac OS.

### EPC 8

Manually and remote controllable patch clamp amplifier.

### EPC 10

The fully computer controlled patch clamp amplifier with built-in LIH 1600 interface.

### PIP 5

Temperature controlled micro pipette puller.

### MIM 4

Motorized 3-axis micromanipulator with video option.

### General notice:

Product names used herein are for identification purposes only and may be trademarks of their respective owners. HEKA disclaims any and all rights in those marks.

We reserve the right to effect technical changes as development progresses. Special versions are available on request. Further technical data are provided by a detailed description, which is available on request. A guarantee of one year applies on all instruments.