

Easy-to-use temperature controller



- Set temperature & GO
- Preset control parameters matched to heater
- Additional input for external sensor

Choose one of these heaters:



Perfusion heaters



Heated microscope stage



Animal heating pads



Objective heater

Why choose npi's temperature controllers:

- Designed for use in sensitive electrophysiological recordings
- **No switching artifacts** due to analog control
- **No electrical noise** due to linear power supply

Contact:

General:
npi electronic GmbH
Phone: +49-7141-9730230
Fax: +49-7141-9730240
sales@npielectronic.com
www.npielectronic.com

North America:
ALA Scientific Instruments
Phone: +1-631-393-6401
Fax: +1-631-393-6407
sales@alascience.com
www.alascience.com